

### Lumipulse® pTau-217/Beta Amyloid 42 Ratio, Plasma New Test

*Notification Date: 10/1/2025*

*Effective Date: 10/7/2025*

The Lumipulse G pTau 217/ $\beta$ -Amyloid 1-42 Plasma Ratio is now available to order. This test is intended to aid healthcare providers to identify patients with amyloid pathology associated with Alzheimer's disease.

The Lumipulse G pTau 217/ $\beta$ -Amyloid 1-42 Plasma Ratio is indicated for adult patients, aged 50 years and older, presenting at a specialized care setting with signs and symptoms of cognitive decline.

Note: Results must be interpreted in conjunction with other patient clinical information. This test is not intended as a screening or stand-alone diagnostic test

#### Ordering Information

**Epic ID:** LAB6535

**Lab Code:** LUMI

#### Specimen Collection Information

**Collection Availability:** Hospital Based Draw Sites and Conrad Jobst (CJT)

**Container:** Lavender-top (EDTA) tube

**Collection Volume:** 3 mL

**Collection Instructions:** Draw full lavender-top (EDTA) tube of blood. Immediately separate plasma from blood via centrifugation (2,000 g for 10 minutes) within two hours after blood collection. Transfer 1mL plasma to Labcorp polypropylene transport tube, seal and **freeze**. Minimum Testing Volume: 0.75 mL plasma  
(Note: This volume does **not** allow for repeat testing.)

#### Stability Requirements

| Temperature             | Period    |
|-------------------------|-----------|
| Room temperature        | 3 hours   |
| Refrigerated            | 1 day     |
| Frozen (-10°C to -60°C) | 1 month   |
| Freeze/thaw cycles      | Stable x0 |

**Storage Instructions:** Frozen

**Causes for Rejection:** Hemolysis; lipemia; heat-treated specimen; gross bacterial contamination; wrong tube type; not stored at proper temperature